




## PRODUCT CATALOGUE

**PARAMOUNT**  
QUEST FOR BETTER HEALTHCARE



## The Paramount Philosophy:

With a vision of being a global leader in the medical disposable industry, Paramount Surgimed believes in providing its customers a timely supply of quality products at a competitive price. The company came into existence in 1993 under the Grover Group of companies capitalizing on the rich and diverse experience of over 60 years. Now, the brand is symbolized with trust and quality in the medical industry. The company started with the manufacturing of Surgical Blades and is now a leading supplier of a huge range of disposable products.

Paramount exports to over 80 countries, has 40 distributors and 400 Dealers across India. Our State-of-the-Art production facility compares to any international standard and is accredited with ISO: 13485:2016, CE Mark and US FDA Device Listing.

## The Paramount Vision:

Paramount Surgimed and all its employees share the vision of their dynamic Managing Director, Mr. Shaily Grover to be a global leader in the Medical Industry. Mr. Grover has ensured that the Group's reputation amongst its associates is synonymous with top quality, updated technology, dedicated approach and customer friendliness.




Corporate Office – Okhla, New Delhi



Factory – Bhiwadi, Rajasthan



Factory – Bawal, Haryana



Factory – Bhiwadi, Rajasthan

Surgical Blades	02
Stitch Cutters	04
Chisel Blades	06
Gouge Blades & Handles	07
Disposable Scalpels	08
Universal Safety Scalpel	09
Biopsy Punch	10
Dermal Curette	12
Ophthalmic Knives	14
Skin Graft Blade & Blood Lancet	16
Disposable Dental Scalpel	17
Myringotomy Knife Blades	18
Face Mask 3 PLY	19
Face Mask N95/FFP2	20
Disposable Adult Diapers	21
Under Pad - AroKleen & PupPee	22
Certificates	23

## Surgical Blades


### Carbon Steel & Stainless Steel

#### Features

- Stainless Steel Surgical Blades and Carbon Steel Surgical Blades
- Packed in peel open pouches with a VCI protective inner line
- Shelf life of 5 years
- Stainless Steel (Grade F as per ISO 7153-1:2016)
- Carbon Steel (Grade- As per BS 2982:1992)
- Complies the Hardness requirement of BS 2982:1992 standard
- Complies the fitment requirement of ISO 7740:1985 standard requirement.



### FITMENT 3



### FITMENT 4







### Details:

Product Name:	Surgical Blades
Size of Inner Box (cm):	8.5 x 8.1 x 2.7
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	47 x 39 x 17
Per 20 feet container:	860 Carton

Size	Carbon Steel	Stainless Steel	Fitment No
9	PCS9	PSS9	3
10	PCS10	PSS10	3
10A	PCS10A	PSS10A	3
11	PCS 11	PSS 11	3
11K	PCS11K	PSS11K	3
11P	PCS11P	PSS11P	3
12	PCS12	PSS12	3
12D	PCS12D	PSS12D	3
13	PCS13	PSS13	3
14	PCS 14	PSS14	3
15	PCS15	PSS15	3
15B	PCS 15B	PSS15B	3
15C	PCS 15C	PSS15C	3
15D	PCS15D	PSS15D	3
15S	PCS15S	PSS15S	3
16	PCS16	PSS16	3
17	PCS17	PSS17	3
18	PCS 18	PSS 18	4
19	PCS 19	PSS 19	4
20	PCS20	PSS20	4
21	PCS 21	PSS 21	4
22	PCS22	PSS22	4
22A	PCS22A	PSS22A	4
23	PCS23	PSS23	4
24	PCS24	PSS24	4
24D	PCS24D	PSS24D	4
25	PCS25	PSS25	4
34	PCS34	PSS34	4
36	PCS36	PSS36	4
36D	PCS36D	PSS36D	4
40	PCS40	PSS40	-
40B	PCS40B	PSS40B	-
60	PCS60	PSS60	8
60B	PCS60B	PSS60B	8

## Post Mortem Blades –Size 60

### Features:

- Used in the fields of Pathology and Histopathology for autopsy
- It is used with handle No.8
- Available in both carbon and stainless steel
- Box of 50
- Individual peel pouch



## Post Mortem Blades –Size 40

### Features:

- The PM40 range of blades are meant to be used with screw type heavy duty handles
- Heavy duty blade
- Box of 10



## Stitch Cutters

### Stainless Steel & Carbon Steel

#### Description

Stitch Cutter is used for general suture removal. Appearing similar to a scalpel, this instrument requires no handle and is a simple and convenient, way to remove basic sutures. The slightly curved, flat, bull nose tip slides under the suture and when rotated, gently removes the suture from the site. Long stitch cutter is made of stainless steel and short stitch cutter is made of carbon steel. The mini stainless steel stitch cutter can be used with handle No.3L, 7 & 9. Alternatively, it is also available on pre-mounted ABS disposable plastic handle.

#### Features

- Five Year Shelf Life – Sterile
- Each blade is protected by VCI paper in a peel-open foil sachet.
- Convenience of single-use disposables.
- Sizes clearly marked on each blade.
- Each box contains 100 pcs

Size	Carbon Steel	Stainless Steel	Handle No.	Handle No.
LSC	–	PLSCSS	PSS3L	PSS9
SSC	PSSCCS	–	PCS3L	PSS9
MSC	–	PMSCSS	PCS3L	PCS7,9

#### Details:

Product Name:	Mini Stitch Cutter
Size of Inner Box (cm):	8.5 x 8.1 x 2.7
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	47 x 39 x 17
Per 20 feet container:	860 Carton



#### Details:

Product Name:	Short Stitch Cutter
Size of Inner Box (cm.):	11 x 8 x 3.5
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm.):	46.5 x 39.8 x 25.5
Per 20 feet container:	570 Carton



**Details:**

Product Name:	Long Stitch Cutter
Size of Inner Box (cm):	14.8 x 7 x 3.1
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	45 x 43 x 34
Per 20 feet container:	430 Carton

**Details:**

Product Name:	Mini Stitch Cutter with Handle
Size of Inner Box (cm):	
Unit/Box:	Bulk Pack of 50 Pcs
Outer Boxes:	Not Using
Export carton:	10 Bags of 50 Pcs x 10 Outer
Export Carton size (cm):	61 x 50 x 3.5
Per 20 feet container:	290 Carton



## Chisel Blades

### Fine Stainless Steel

#### Description

FINE (CHISEL) Surgical Blades Stainless Steel are available in the following sizes: Nos. P-SM61, P-SM61S, P-SM62, P-SM62S, P-SM63, P-SM65, P-SM65A, P-SM67, P-SM681, P-SM69, P-SM90, P-SM91 Each blade is protected by VCI paper in a peel-open foil sachet. Sterilized by Gamma Radiation of 2.5 M. Rads.

#### Details:

Product Name:	Chisel Blades (Fine)
Size of Inner Box (cm):	82 x 32 x 22
Unit/Box:	25 Pcs
Outer Boxes:	25 Boxes x 25 Units
Export carton:	25 Boxes x 10 Outer
Export Carton size (cm):	47 x 39 x 17
Per 20 feet container:	860 Carton

Size	Product Code	Handle
61	PSM61	PSMH
61V	PSM61V	PSMH
62	PSM62	PSMH
62V	PSM62V	PSMH
63	PSM63	PSMH
64	PSM64	PSMH
65	PSM65	PSMH
66	PSM66	PSMH
67	PSM67	PSMH
68	PSM68	PSMH
69	PSM69	PSMH
69B	PSM69B	PSMH
90	PSM90	PSMH
91	PSM91	PSMH






## Gouge Blades & Handles

### Carbon Steel

#### Description

These surgical blades are used for feet surgeries and podiatry purposes i.e. manicure and pedicure. Gouge Blades are used by various podologists for hand, feet and skin care. They are available in the following sizes: Nos. 1, 2, 3, 4, 5, 6, 8, 10, 12, 15 and 18. Each blade is protected by VCI paper in a peel- open foil sachet. Gouge blades fit into handle no. 3 and pedal gouge handles.



#### Details:

Product Name:	Gouge Blades No. 1 to 5
Size of Inner Box (cm):	10 x 10 x 3.3
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	59 x 45 x 21
Per 20 feet container:	490 Carton

Size	Product Code	Size	Product Code	Fitment No	Handle
1	PGCS1	1V	PGCS1V	3	PGMH
2	PGCS2	2V	PGCS2V	3	PGMH
3	PGCS3	3V	PGCS3V	3	PGMH
4	PGCS4	1R	PGCS1R	3	PGMH
5	PGCS5	2R	PGCS2R	3	PGMH
6	PGCS6	3R	PGCS3R	3	PGMH
8	PGCS8			3	PGMH
10	PGCS10			3	PGMH
12	PGCS12			3	PGMH
15	PGCS15			3	PGMH
18	PGCS18			3	PGMH

#### Details:

Product Name:	Gouge Blades No. 6 to 18
Size of Inner Box (cm):	13.8 x 10 x 3.3
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	58 x 40 x 29.5
Per 20 feet container:	400 Carton




## Disposable Scalpels

### Stainless Steel blade

#### Features


- Pre-assembled stainless steel surgical blades mounted on ABS plastic handle.
- The device is used to prevent injury during fixing the blade on the metal handle.
- Fluted finger grips and the blades are safely protected by LDPE guard.
- Individually packaged in a peel pouch.
- Convenience of single-use.
- Size clearly marked on each scalpel blade.
- Each dispenser box contains 10 scalpels.
- Shelf life of 5 years



#### Details:

Product Name:	Disposable Scalpel
Size of Inner Box (cm):	20 x 6.8 x 3.5
Unit/Box:	10 Pcs
Outer Boxes:	10 Boxes x 10 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	44.5 x 42.5 x 38
Per 20 feet container:	370 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PDSS9	3	18	PDSS18	4
10	PDSS10	3	19	PDSS19	4
10A	PDSS10A	3	20	PDSS20	4
11	PDSS11	3	21	PDSS21	4
11K	PDSS11K	3	22	PDSS22	4
12	PDSS12	3	22A	PDSS22A	4
13	PDSS13	3	23	PDSS23	4
14	PDSS14	3	24	PDSS24	4
15	PDSS15	3	25	PDSS25	4
15B	PDSS15B	3	36	PDSS36	4
15C	PDSS15C	3			
16	PDSS16	3			
17	PDSS17	3			



## Universal Safety Scalpel

### Stainless Steel blade

#### Features

- **Safety Mode:** When you retract back after fully extending the blade, it makes a click sound to denote it is in safety mode. Safety mode means that the scalpel is not fully extended. This feature avoids injuries when the blade is half way extended as it retracts back and does not pierce through. It has to be fully extended to use.
- **Safety Lock:** if you completely retract back it makes a loud click denoting it is in safety lock position. Safety lock position means that the blade cannot be reused once it is in lock position. This feature forces the doctor to dispose the product and minimizes the risk of transmissibility of blood-borne infectious agents.
- Compatible for a left handed as well as right handed surgeon.
- Safety feature in pre use and post use of product.
- One-handed activation allowing the clinician to focus on the surgical site.
- The retractable design protects the blade from exposure, minimizing the risk of accidental cuts.
- Universal size handle that fits all types of hand sizes, and ideal for kits, trays & packs.
- One handle size fits all blade sizes.

#### Details:

Product Name:	Safety Scalpel
Size of Inner Box (cm):	10 Pcs
Unit/Box:	17 x 9 x 4
Outer Boxes:	5 Boxes x 10 Units
Export carton:	5 Boxes x 20 Outer
Export Carton size (cm):	54 x 39 x 44.5
Per 20 feet container:	280 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PSF9	3	18	PSF18	4
10	PSF10	3	19	PSF19	4
10A	PSF10A	3	20	PSF20	4
11	PSF11	3	21	PSF21	4
11K	PSF11K	3	22	PSF22	4
13	PSF13	3	22A	PSF22A	4
14	PSF14	3	23	PSF 23	4
15	PSF15	3	24	PSF24	4
15B	PSF15B	3	25	PSF25	4
15C	PSF15C	3			
16	PSF16	3			
17	PSF17	3			



## Biopsy Punch

### Stainless Steel blade

#### Description

Biopsy punch is a sharp rounded blade tip used to remove tissue to be tested during procedures. The instrument is sterile with a round plastic handle. The punch is available in a variety of tip sizes and is sold in both sterile single pack as well as in bulk quantities.


Finely ground seamless, stainless steel blade ensures safe, clean and easy sample extraction. The sleek ribbed handle clear view of the biopsy site, comfort and control.

Virgin grade ABS plastic provides strength and sturdiness.

#### Details:

Product Name:	Biopsy Punch
Size of Inner Box (cm):	17 x 9 x 3.5
Unit/Box:	10 Pcs
Outer Boxes:	10 Boxes x 10 Units
Export carton:	10 Boxes x 5 Outer
Export Carton size (cm):	40 x 20 x 51.5
Per 20 feet container:	660 Carton


Size mm	Product Code	Product Code	Product Code
1.0	BPS1.0	BPK1.0	BPP1.0
1.5	BPS1.5	BPK1.5	BPP1.5
2.0	BPS2.0	BPK2.0	BPP2.0
2.5	BPS2.5	BPK2.5	BPP2.5
3.0	BPS3.0	BPK3.0	BPP3.0
3.5	BPS3.5	BPK3.5	BPP3.5
4.0	BPS4.0	BPK4.0	BPP4.0
5.0	BPS5.0	BPK5.0	BPP5.0
6.0	BPS6.0	BPK6.0	BPP6.0
7.0	BPS7.0	BPK7.0	BPP7.0
8.0	BPS8.0	BPK8.0	BPP8.0
10.0	BPS10.0	BPK10.0	BPP10.0
12.0	BPS12.0	BPK12.0	BPP12.0
15.0	BPS15.0	BPK15.0	BPP15.0





### BPS Features

- Ergonomically designed for precise skin excisions
- Unique, ribbed, easy-grip barrel and pressure assistance platform
- Size embossed on each punch for easy identification.
- Sterile and ready to use.
- Wide range of sizes



### BPK Features

- Seamless stainless steel sharp cutting edge.
- Protective cap for each punch to maintain sharpness
- Hollow handle provides easy specimen retrieval
- Minimal damage to tissues
- Wide range of sizes



### BPP Features

- Premium quality punch
- Hard blister packaging
- Option to print custom brand on handle
- Pressure assistance platform and ribbed handle
- Micro grinding finish
- Size embossed on handle



## Dermal Curette

### Stainless Steel blade

#### Description

The curette has a firm place in the range of dermatological, diagnostic and therapeutic instruments. It is one of the most versatile instruments. For that reason, it is increasingly also used in other disciplines, such as plastic surgery and gynaecology.


The aim of a curette is the removal of a superficial skin changes with minimal damage to the surrounding tissue.

The curette is extremely sharp thus ensuring precise and gentle tissue removal. It can also be used for the treatment of superficial malignant tumours. The outstanding cutting ability makes working with it much easier compared with conventional products.

#### Details:

Product Name:	Dermal Curette
Size of Inner Box (cm):	21 x 7 x 4
Unit/Box:	10 Pcs
Outer Boxes:	5 Boxes x 10 Units
Export carton:	5 Boxes x 10 Outer
Export Carton size (cm):	46.5 x 46 x 25
Per 20 feet container:	500 Carton

Size mm	Product Code	Product Code	Product Code	Product Code
2.0	DCS2.0	DCP2.0	DCK2.0	DDC2.0
3.0	DCS3.0	DCP3.0	DCK3.0	DDC3.0
4.0	DCS4.0	DCP4.0	DCK4.0	DDC4.0
5.0	DCS5.0	DCP5.0	DCK5.0	DDC5.0
7.0	DCS7.0	DCP7.0	DCK7.0	DDC7.0




### DCS Features

- Seamless stainless steel sharp cutting edge
- Easy dynamic gripping handle
- Sterile and ready to use
- Sharp side indicator
- Wide range of sizes




### DCK Features

- Sturdy & Durable Stainless Steel
- Symmetrically Ground Beveled Edge
- Ultra Sharp Inside Grinding
- Sharp side indicator
- Ribbed on all sides
- Triangular Handle
- Sterile and ready to use
- Wide range of sizes




### DCP Features

- Long 'Pen Type' Grip
- Branding Labeling area
- Ribbed handle for easy control
- Minimum damage to surrounding tissues
- Sharp side indicator
- Sterile and ready to use
- Wide range of sizes



### DDC Features

- Seamless stainless steel sharp cutting edge
- One in metal with a sharp edge
- Sharp side indicator
- Covered by a protective cap
- Ultra Sharp Inside Grinding
- Sterile and ready to use
- Wide range of sizes



## Ophthalmic Knives

### Round-Wire Blades



### Slit-Blades (Keratome) For Phaco Stab Incision

#### Product No. Description

P-912301 2.3 MM Sharp Tip Straight  
P-912501 2.5 MM Sharp Tip Straight  
P-912601 2.65 MM Sharp Tip Straight  
P-912801 2.8 MM Sharp Tip Straight  
P-912901 3.0 MM Sharp Tip Straight  
P-913201 3.2 MM Sharp Tip Straight  
P-913501 3.5 MM Sharp Tip Straight  
P-915101 5.1 MM Sharp Tip Straight

#### Product No. Description

P-912361 2.3 MM Sharp Tip Angled 45° 19 Gauge  
P-912561 2.5 MM Sharp Tip Angled 45° 19 Gauge  
P-912661 2.65 MM Sharp Tip Angled 45° 19 Gauge  
P-912861 2.8 MM Sharp Tip Angled 45° 19 Gauge  
P-912961 3.0 MM Sharp Tip Angled 45° 19 Gauge  
P-913261 3.2 MM Sharp Tip Angled 45° 19 Gauge  
P-913561 3.5 MM Sharp Tip Angled 45° 19 Gauge  
P-915061 5.1 MM Sharp Tip Angled 45° 19 Gauge

#### Double Bevel

P-912808 2.8 MM Sharp Tip Straight  
P-912908 3.0 MM Sharp Tip Straight  
P-913208 3.2 MM Sharp Tip Straight

#### Double Bevel

P-912868 2.8 MM Sharp Tip Angled 45° 19 Gauge  
P-912968 3.0 MM Sharp Tip Angled 45° 19 Gauge  
P-913268 3.2 MM Sharp Tip Angled 45° 19 Gauge



### Slit-Blades (Keratome) For I.O.L. Enlargers

#### Product No. Description

P-914001 4.0 MM Sharp Tip Straight  
P-915201 5.2 MM Sharp Tip Straight  
P-915501 5.5 MM Sharp Tip Straight  
P-916001 6.0 MM Sharp Tip Straight  
P-916201 6.2 MM Sharp Tip Straight

#### Product No. Description

P-914061 4.0 MM Sharp Tip Angled 45° 16 Gauge  
P-915561 5.5 MM Sharp Tip Angled 45° 15 Gauge  
P-916061 6.0 MM Sharp Tip Angled 45° 15 Gauge  
P-916261 6.2 MM Sharp Tip Angled 45° 15 Gauge



### Crescent Blades For Groove Or Pocket Or Tunnel Incision

#### Product No. Description

P-950001 Crescent Sharp Tip Straight

#### Product No. Description

P-950002 Crescent Bevel Up Angled 45° 20 Gauge  
P-950003 Crescent Bevel Down Angled 45° 20 Gauge  
P-950004 Crescent Angled 45° 20 Gauge

#### Double Bevel

P-950005 Sharp Tip Straight



### Lance Tip Blades For Initial Incision

#### Product No. Description

P-931501 15° Degree 24 Gauge  
P-933001 30° Degree 24 Gauge  
P-934501 45° Degree 24 Gauge

#### Product No. Description

P-913501 15° Degree 24 Gauge  
P-915101 30° Degree 24 Gauge  
P-915161 45° Degree 24 Gauge

#### Product No. Description

P-914101 15° Degree 24 Gauge  
P-914161 30° Degree 24 Gauge  
P-915261 45° Degree 24 Gauge



## Round-Wire Blades



## MVR Blades For Water Tight Self Sealing Incision

### Product No. Description

P-975559 MVR 24 Gauge

P-975560 MVR 20 Gauge

P-975561 MVR 19 Gauge

### Product No. Description

AC STILETIO BLADES

For Blumenthal or Side Port Incision

P-985560 MVR 20 Gauge Angled 45°

P-985561 MVR 19 Gauge Angled 45°



## Miniature Blades Alternative To General Surgery For Precise Incision

### Product No. Description

SPOON BLADES

For Groove or Pocket or Tunnel Incision

### Product No. Description

P-6820 1.5 MM Angled 45°

P-6821 3.0 MM Angled 60°

P-6821E 3.0 MM Extended Neck Angled 60°



## Miniature Blades Alternative To General Surgery For Precise Incision


### Product No. Description

SCLERAL BLADES

For Limbal and Scleral Incision

P-5700 Wide Tip Blade

P-6821 Narrow Tip Blade



### Details:

Product Name:	Ophthalmic Knives
Size of Inner Box (cm):	20 x 5 x 4
Unit/Box:	6 Pcs
Outer Boxes:	25 Boxes x 6 Units
Export carton:	25 Boxes x 2 Outer
Export Carton size (cm):	44 x 28.5 x 23
Per 20 feet container:	1000 Carton

## Skin Graft Blade

### Carbon Steel

#### Description

##### Skin Graft Blade, Simplex, Duplex


Surgical Skin Graft Blades are 158 mm long, individually wrapped in VCI paper in single peel packets in packs of 10.

These blades can be used with Watson, Braithwaite or Cobbett Knives.

#### Details:

Product Name:	Skin Graft Blades
Size of Inner Box (cm):	20.7 x 3.6 x 2
Unit/Box:	10 Pcs
Outer Boxes:	10 Boxes x 10 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	44 x 23 x 22.5
Per 20 feet container:	1250 Carton

Size	Carbon Steel	Stainless Steel
SG	PCSG	PSSG
SGS	PCSGS	PSSGD
SGD	PCSGS	PSSGD



## Blood Lancet

#### Description

Supplied sterile and individually packed. Lancet comes in a pack of 200 pieces per box. the tip is specially designed for easy pricking, which minimizes patient trauma.

#### Details:

Product Name:	Blood Lancet
Product Code:	PSSBL
Size of Inner Box (cm.):	12.5 x 7 x 2.5
Unit/Box:	200
Export carton:	100 Boxes x 200 Units
Export Carton size (cm.):	38.5 x 26.5 x 15.5
Per 20 feet container:	1500 Cartons



## Disposable Dental Scalpel

### Stainless Steel

#### Description

Dental disposable scalpels initially conceived for parodontology and Stomatology, because their form allows to do precise and fine Incision, even in hardly accessible areas.


These instruments have numerous applications in other disciplines like cardiology, plastic surgery, gynecology and neurology.

They can be presented in two ways: with a plastic molded handle, or separately, in order to be fixed on a special handle.

#### Details:

Product Name:	Dental Scalpel
Size of Inner Box (cm):	23 x 7 x 5.6
Unit/Box:	5 Pcs
Outer Boxes:	5 Boxes x 5 Units
Export carton:	5 Boxes x 10 Outer
Export Carton size (cm):	61 x 50 x 3.5
Per 20 feet container:	290 Carton

Size	Product Code
1	DBSS1
2	DBSS2
3	DBSS3



## Myringotomy Knife Blades

### Description

Myringotomy is a surgical procedure of the eardrum or tympanic membrane. The procedure is performed by making a small incision with a myringotomy knife through the layers of tympanic membrane


### Features:

- Optimised micrometre-accurate precision grinding
- Double edged allowing for cutting movements in both directions
- Ultra-sharp cutting edges of superb durability
- Consistently sharp, angled blade facilitates a clean cut
- Spear tip, Narrow shaft
- Lance tip, Downward cut, Narrow shaft

### Details:



Product Name:	Myringotomy Knife Blade
Size of Inner Box (cm):	24.2 x 8 x 2.5
Unit/Box:	5 Pcs
Outer Boxes:	--
Export carton:	--
Export Carton size (cm):	--
Per 20 feet container:	--

Size	Product Code
Lance	PSSML
Spear	PSSMS





## Face Mask - 3 Ply



### 3 PLY Mask - Non Surgical



Smooth Fiber  
Material



SMS Fabric

### 3 PLY Mask - Surgical



Meltblown  
Fabric



UV Cleaned

## FEATURES



Adjustable Nose  
Clip



Ultrasonic  
Sealing



One Size  
Fits Most



Seamless  
Crimping



Comfortable  
& Breathable



Quality  
Material



3 Layer  
Protection



Anti Droplets

### Details:

Product Name:	3 PLY Mask - Non Surgical
Product Code	P03N
Size of Inner Box (cm.):	23 x 11 x 19.5
Unit/Box:	100
Export Carton size (cm.):	57 x 48 x 41
Per 20 feet container:	230 Carton

### Details:

Product Name:	3 PLY Mask - Surgical
Product Code	P03M
Size of Inner Box (cm.):	23 x 10.5 x 10
Unit/Box:	50
Export Carton size (cm.):	55 x 48 x 42
Per 20 feet container:	230 Carton

## Face Mask - N95/FFP2



without Respirator



with Respirator



UV Cleaned



Filters Particulates



BFE ≥ 95%



PFE ≥ 95%

### FEATURES



Adjustable Noise Clip



Soft Elastic Ear Loops



High Quality Breathable  
Material and Valve



Seamless Flow of  
Exhaled Air

#### 1st Layer

Fine mesh.  
Filters unwanted particulates

#### 3rd Layer

Hot air mesh.  
Blocks smaller particulates

#### 2nd Layer

Heavy meltblown filter mesh.  
Blocks odour, virus, pollution

#### 4th Layer

Extra meltblown.  
Filters remaining particulates

#### 5th Layer

Extra soft non woven  
layer to fit the face

### Details:

Product Name:	N95/FFP2 without Respirator
Product Code:	P95W
Size of Inner Box (cm.):	30.5 x 13.5 x 17.5
Unit/Box:	50
Export Carton size (cm.):	63.5 x 56 x 55
Per 20 feet container:	140 Carton

### Details:

Product Name:	N95/FFP2 with Respirator
Product Code:	P95R
Size of Inner Box (cm.):	24.5 x 13.5 x 17.5
Unit/Box:	25
Export Carton size (cm.):	51 x 42 x 54.5
Per 20 feet container:	220 Carton

## Disposable Adult Diapers


### For Urinary Incontinence Management

#### Features:

- Dry-Fast Core, Super absorbent.
- Moisture repellent.
- Security barriers.
- Waist Elastication and curved leg elastic.
- Superfit Band Wetness Indicator.
- Three absorbency levels and three sizes.
- Compact and convenient packaging.
- Also available in options of double cotton core & ADL.

#### Benefits:

- Faster absorption, drier surface.
- Better containment of faces.
- Reduces urine leakage.
- Improves comfort.
- Comfortable fit.
- Easier to fix and adjust tapes.
- A Guide for carers.
- Meets individual requirements with cost effective care.
- Easy to handle less storage space required.



#### Details:

Product Name:	Adult Diapers (Medium)
Size of Poly Bag (cm.):	20 x 23.0 x 17.5
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	46 x 35.5 x 60
Per 20 feet container:	275 Cartons

#### Details:

Product Name:	Adult Diapers (Large)
Size of Poly Bag (cm.):	22 x 24 x 17.5
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	48.5 x 35.5 x 66
Per 20 feet container:	250 Cartons

#### Details:

Product Name:	Adult Diapers (Extra Large)
Size of Poly Bag (cm.):	23.5 x 24 x 17.5
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	49 x 35.5 x 70
Per 20 feet container:	220 Cartons



## Under Pads

### AroKleen

#### Description

- **Anti-Leakage Design:**  
The product has extra width and length designed for better absorption.
- **Dry And Comfortable:**  
It has large absorbing padding using the reticular embossed design to introduce the water diversion towards every side keep the surface comfortable.
- **Skin Friendly:**  
Using the super-nonwoven as the top layer, it feels comfortable and prevents skin rashes.
- **International Quality:**  
Strict controls of product material and production assuring safety and health.


#### Details:

Product Name:	Under Pad
Size of Poly Bag (cm.):	29 x 20 x 10
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	41 x 31 x 59
Per 20 feet container:	350 Cartons



### PupPee –Pet Training Pads

#### 4 LAYERED LEAK PROOF TECHNOLOGY!



#### The Advantages of Using PupPee Pads

Our pads provide worry-free training for your puppy & peace of mind for you. They also can be used as an alternative to outdoor "relief" for your adult dog.

- Locks in wetness with quilted top layer
- Ensures maximum absorption power with 2 moisture-locking layers
- Protects floor & carpets with leak-proof backing
- Provides easy training with built-in attractant
- Protect your floors
- Easy cleanup & disposal
- Quilted for extra absorbency
- Less accidents means less disciplining
- Dog attractant reinforces training
- Durable, tear resistant



**DNV·GL**

## EC CERTIFICATE

### Full Quality Assurance System

Certificate No.: 10000400574-PA-NA-IND      Project No.: PRIC-213992-2010-PRC-IND      Valid Until: 27-05-2024

This is to certify that the quality system of:

**Paramount Surgimed Ltd.**  
Works: A-106, RIICO Industrial Area, Bhiwadi – 301 019, District Alwar, Rajasthan, India

For design, production and final product inspection/testing of:  
**STERILE SURGICAL DISPOSABLE MEDICAL DEVICES**

Has been assessed with respect to:  
**THE CONFORMITY ASSESSMENT PROCEDURE DESCRIBED IN ANNEX II EXCLUDING SECTION 4 OF COUNCIL DIRECTIVE 93/42/EEC ON MEDICAL DEVICES, AS AMENDED**

and found to comply.

Further details of the product(s) and conditions for certification are given overleaf.

Place and date:  
Havik, 15 September 2020

For:  
**DNV GL PRESAFE AS**  
Notified Body No.: 2460

Mariann Jeremiassen

The certificate is digitally verified by blockchain technology. For more info, visit [www.dnv-gl.com/resources/certificates-on-the-blockchain.html](http://www.dnv-gl.com/resources/certificates-on-the-blockchain.html)

NOTIFIED BODY: DNV GL PRESAFE AS, Vertovæien 3, N-1363 Havik, Norway - Registered Enterprise No. NO 987 987 402 MVA

ISSD C0-078-A Rev 0.0      Page 1 of 4

EC Certificate

**DNV·GL**

Certificate No.: 10000400574-PA-NA-IND      Project No.: PRIC-213992-2010-PRC-IND      Valid Until: 27-05-2024

**Jurisdiction**

Application of Council Directive 93/42/EEC of 14 June 1993, adopted as "Forskrift om Medisinsk Utstyr" by the Norwegian Ministry of Health and Care Services.

**Certificate history:**

Revision	Description	Issue Date
0.0	Original Certificate	2020-09-15

Products covered by this Certificate:

Product Description	Product Name	Class
Sterile Curette Dermal	Curette Dermal Size: 2, 3, 4, 5 & 7	Ila
Sterile Stich cutters in Carbon Steel and Stainless Steel	Stitch Cutters Long, Short, Mini	Ila
Sterile Surgical Blades in Carbon Steel and Stainless Steel	Surgical blades 1, 2, 3, 4, 5, 6, 8, 9, 10, 10A, 11, 11K, 12, 13, 14, 15, 15B, 15C, 15D, 15S, 16, 17, 18, 19, 20, 21, 22, 22A, 23, 24, 25, 36, 40, 40B, 60, 60B, 1R, 2R, 3R, 1V, 2V, 3V, 11P, 12D, 24D, 34, 36D	Ila
Sterile Disposable Scalpels in Carbon Steel and Stainless Steel	Disposable Scalpels 1, 2, 3, 4, 5, 6, 8, 9, 10, 10A, 11, 11K, 12, 13, 14, 15, 15B, 15C, 15D, 15S, 16, 17, 18, 19, 20, 21, 22, 22A, 23, 24, 25, 36, 1R, 2R, 3R, 1V, 2V, 3V, 11P, 12D, 24D, 34, 36D	Ila
Sterile Safety Scalpels in Carbon Steel and Stainless Steel	Safety Scalpels 6, 8, 10, 10A, 11, 11K, 12, 13, 14, 15, 15B, 15C, 15D, 15S, 16, 17, 18, 19, 20, 21, 22, 22A, 23, 24, 25, 36, 1R, 2R, 3R, 1V, 2V, 3V, 11P, 12D, 24D, 34, 36D	Ila
Sterile Fine Blades / Chisel Blade / Microsurgery blades	61, 62, 63, 64, 65, 66, 67, 68, 69, 69B, 61V, 62V, 90, 91	Ila
Sterile Ophthalmic Blades	Keratome: P-912301, P-912501, P-912601, P-912801, P-912901, P-913201, P-913501, P-915001, P-912361, P-912561, P-912661, P-912861, P-912961, P-913261, P-913561, P-915061, P-912808, P-912908, P-913208, P-912868, P-912968, P-913268, P-914001, P-915201	Ila

NOTIFIED BODY: DNV GL PRESAFE AS, Vertovæien 3, N-1363 Havik, Norway - Registered Enterprise No. NO 987 987 402 MVA

ISSD C0-078-A Rev 0.0      Page 2 of 4

EC Certificate

**DNV·GL**

Certificate No.: 10000400574-PA-NA-IND      Project No.: PRIC-213992-2010-PRC-IND      Valid Until: 27-05-2024

P915501, P-916001, P-916201, P-914061, P-915161, P915561, P-916061, P-916261 Crescent: P-950001, P-950002, P-950003, P-950004, P950005 Lance Tip: P-931501, P-933001, P-934501, P-913501, P915101, P-915161, P-914101, P-914161, P-915261 MVR: P-975559, P-975560, P-975561, P-985560, P-985561 Spoon: P-6820, P-6821, P-6821E Scleral: P-5700, P5710		
Sterile Biopsy Punches	Biopsy Punches: 1mm, 1.5mm, 2mm, 2.5mm, 3mm, 3.5mm, 4mm, 5mm, 6mm, 7mm, 8mm, 10mm, 12mm, 15mm	Ila
Sterile Skin Graft Blades	Simplex, Duplex	Ila
Sterile Blood Lancet	Standard	Ila
Sterile Myringotomy	Lance and Spear	Ila

The complete list of devices is filed with the Notified Body

**Sites covered by this certificate**

Site Name	Address
Paramount Surgimed Ltd.	A-106, RIICO Industrial Area, Bhiwadi – 301 019, District Alwar, Rajasthan, India

**EU Representative**  
Medical Device Safety Service GmbH,  
Schiffgraben 41, D-30175, Hannover, Germany.

NOTIFIED BODY: DNV GL PRESAFE AS, Vertovæien 3, N-1363 Havik, Norway - Registered Enterprise No. NO 987 987 402 MVA

ISSD C0-078-A Rev 0.0      Page 3 of 4

EC Certificate

**DNV·GL**

Certificate No.: 10000400574-PA-NA-IND      Project No.: PRIC-213992-2010-PRC-IND      Valid Until: 27-05-2024

**Terms and conditions**

The certificate is subject to the following terms and conditions:

- Any producer (see 2001/95/EC for a precise definition) is liable for damage caused by a defect in his product(s), in accordance with directive 85/374/EEC, as amended, concerning liability of defective products.
- The certificate is only valid for the products and/or manufacturing premises listed above.
- The Manufacturer shall fulfil the obligations arising out of the quality system as approved and uphold it so that it remains adequate and efficient.
- The Manufacturer shall inform Presafe of any intended updating of the quality system and Presafe will assess the changes and decide if the certificate remains valid.
- Periodical audits will be held, in order to verify that the Manufacturer maintains and applies the quality system. Presafe reserves the right, on a spot basis or based on suspicion, to pay unannounced visits.

The following may render this Certificate invalid:

- Changes in the quality system affecting production.
- Periodical audits not held within the allowed time window.

**Conformity declaration and marking of product**

When meeting with the terms and conditions above, the producer may draw up an EC declaration of conformity and legally affix the CE mark followed by the Notified Body identification number of Presafe.

End of Certificate

NOTIFIED BODY: DNV GL PRESAFE AS, Vertovæien 3, N-1363 Havik, Norway - Registered Enterprise No. NO 987 987 402 MVA

ISSD C0-078-A Rev 0.0      Page 4 of 4

EC Certificate






ISO 9001 : 2015



ISO 14001 : 2015



ISO : 13485



ISO : 13485

# CERTIFICATE OF REGISTRATION

This is to certify that the management system of:

## Paramount Surgimed Limited

Main Site: 1, L.S.C., Okhla Industrial Area, Phase-II, New Delhi, 110020, India (FIN F004065)

Additional Site1: A-106, RIICO Industrial Area, Bhiwadi, Alwar, 301019, Rajasthan, India (FIN F004068)

Additional Site2: Plot No. 33, Sector-14, Phase-2, IMT Bawal, 123501, Haryana, India (FIN F004067)

has been registered by Intertek, an MDSAP recognized auditing organization, as conforming to the requirements of:

### ISO 13485:2016

**Australia:** Therapeutic Goods (Medical Devices) Regulations, 2002, Schedule 3 Part 1 (excluding Part 1.6)

**Brazil:** Federal Law n. 6360/76; RDC ANVISA n. 16/2013; RDC ANVISA n. 23/2012; RDC ANVISA n. 67/2009; RDC ANVISA n. 56/2001

**Canada:** Medical Devices Regulations – Part 1- SOR 98/282

**United States:** 21 CFR 820, 21 CFR 803, 21 CFR 806, 21 CFR 807 (Subparts A to D)

### The management system is applicable to:

*Design and Development, Manufacture and Supply of Surgical Blades, Disposable Scalpels & Disposable Safety Scalpels (Sterile & Non Sterile in Carbon Steel and Stainless Steel), Blood Lancets, Ophthalmic Blades, Skin Graft Blades, Curette Dermal and Biopsy Punches, Fine Blades/Chisel Blade/Microsurgery blades and Myringotomy sterile.*

*Head Office: Sales and Marketing, Management, Purchasing*

*Additional Site1: Design and Development, Manufacturing, Process Control, Quality Assurance.*

*Additional Site2: Manufacturing of Plastic components and process control.*

Certificate Number:

0097007

Initial Certification Date:

2019-12-17

Certification Effective Date:


2019-12-17

Certification Expiry Date:

2022-12-16



Intertek

  
Calin Moldoveanu

President, Business Assurance

Intertek Testing Services NA, Inc.  
900 Chelmsford Street  
Lowell, MA, USA 01851



In the issuance of this certificate, Intertek assumes no liability to any party other than to the Client, and then only in accordance with the agreed upon Certification Agreement. This certificate's validity is subject to the organization maintaining their system in accordance with Intertek's requirements for systems certification. Validity may be confirmed via email at [certificate.validation@intertek.com](mailto:certificate.validation@intertek.com) or by scanning the code to the right with a smartphone. The certificate remains the property of Intertek, to whom it must be returned upon request. Validity of this certificate may be verified at <http://www.intertek.com/business-assurance/certificate-validation/>

CT-MDSAP-2016-NA-EN-LT-P-30.apr.18



MDSAP



1, L.S.C., Okhla Industrial Area, Phase-II, New Delhi-110 020, India

Tel +91-11-46436600, 46436601, 46436666, 26389812, 13

Fax +91-11-26389815, 41616555

[sales@paramountblades.com](mailto:sales@paramountblades.com)

[export2@paramountblades.com](mailto:export2@paramountblades.com)

[www.paramountsurgimed.com](http://www.paramountsurgimed.com)

(Disclaimer- All measurements in weight & size are indicative for reference only.)

PSL/CTLG/01, Issue No. 01, Issue Date : 26/11/2020, Rev. No. 00, Rev. Date : NIL,